

How to connect a chatwrap case

This video will show you how to connect your iPad to a chatwrap case via Bluetooth. The first thing you will need to do is make sure that both your iPad and the chatwrap case have enough charge. If the case is not charged (or not charged enough) it will not appear in available devices for connection options. After everything is charged, go to your iPad's settings. From there, find the Bluetooth option from the menu on the left side of the screen. Tap to turn the Bluetooth "ON". From the drop down menu of available devices choose the one that says "Chatwrap". After a few moments a message will appear next to the speechcase's name saying „connected“.

You can also test if it worked by opening grid and "speak" a sentence to see where if the sound is coming from the iPad's or the case's speaker.