

24 HOUR URINE COLLECTION

To collect a 24 hour urine sample correctly, the patient will need:

- 1) **The appropriate coloured 24 hour urine container**
available from www.nbt.nhs.uk/pathologyconsumables



Plain

Albumin
Bence Protein
Cortisol
Copper
Electrolytes
Protein

HCl Acid

Calcium (Ca)
Citrate
Magnesium
Metadrenalines
Oxalate (Ox)
Phosphate
CaOx stone former**
CaOx stone former**
(Follow Up)
Cystine stone former**

Acetic Acid

5 HIAA

Alkaline

Urate
UA stone former

Acid Day & Night

Cystine**
(Day and Night)

**** Hospital Requests Only**

- 2) Information to complete the sample label, including:
 - a. An **ICE request label**,
 - b. Their **NHS number**
 - c. **Name of the test** required

- 3) **Instruction leaflets**, which should arrive with the container, are also available from www.nbt.nhs.uk/severn-pathology/patient-information

**IF DETAILS ARE NOT FULLY COMPLETED ON THE SAMPLE,
SAMPLES WILL NOT BE PROCESSED**