

24 HOUR URINE COLLECTION

To collect a 24 hour urine sample correctly, the patient will need:

- 1) **The appropriate coloured 24 hour urine container** available from www.nbt.nhs.uk/pathologyconsumables



Plain	HCl Acid	Acetic Acid	Alkaline	Acid Day & Night
Albumin Bence Protein Cortisol Copper Electrolyte Protein	Calcium (Ca) Citrate Magnesium Metadrenalines Oxalate (Ox) Phosphate CaOx stone former** CaOx stone former** (follow-up) Cystine stone former**	5 HIAA	Urate UA stone former	Cystine** (Day and Night)

**** Hospital Requests Only**

- 2) Information to complete the sample label, including:
 - a. An **ICE request label**,
 - b. Their **NHS number**
 - c. **Name of the test** required

- 3) **Instruction leaflets**, which should arrive with the container, are also available from www.nbt.nhs.uk/severn-pathology/patient-information

**IF DETAILS ARE NOT FULLY COMPLETED ON THE SAMPLE,
SAMPLES WILL NOT BE PROCESSED**